

DUAL COMPARTMENT WARMING CABINET

Introduction

SKYTRON J3 Series warming cabinets are designed to maintain clinically ideal heated storage of blankets, sterile water, and saline solutions.

The body of the warming cabinet is double walled stainless steel construction with insulation for increased heating efficiency. Doors are made of double walled stainless steel (J3 models) or double paned tempered glass framed with extruded aluminum (J3G models).

The SS2207-J3(G) has a total interior capacity of 17.3 cubic feet, from a top chamber measuring 5.1 cubic feet and a lower chamber measuring 12.2 cubic feet.

Standard Features

- Temperature Range 90°-160°F (32°-71°C)
- Digital Thermostatic Control and Display
- Lock-out Capability (Temperature Setting)
- Spring Loaded Hinged Doors
- Front Access Serviceability
- Trend Analysis w/ USB Retrieval


Options

- Solid Door Configuration (SS2207-J3)
- Glass Door Configuration (SS2207-J3G)
- Caster Kit (F2-012-90)



Glass Door Dual Compartment Warming Cabinet
(SS2207-J3G)

Technical Drawing



Technical Data

Model:	SS2207-J3	SS2207-J3(G)
Exterior Dimensions		
Width:	30" (762mm)	30" (762mm)
Height:	74" (1880mm)	74" (1880mm)
Depth:	24" (610mm)	24" (610mm)
Interior Dimensions Upper Chamber		
Width:	24 7/8" (632mm)	24 7/8" (632mm)
Height:	15 1/4" (387mm)	15 1/4" (387mm)
Depth:	23 3/8" (594mm)	23 3/8" (594mm)
Interior Dimensions Lower Chamber		
Width:	24 7/8" (632mm)	24 7/8" (632mm)
Height:	36 3/8" (924mm)	36 3/8" (924mm)
Depth:	23 3/8" (594mm)	23 3/8" (594mm)
Electrical		
	120 VAC, 50/60 Hz, 1 PH, 10.8 A	120 VAC, 50/60 Hz, 1 PH, 10.8 A
Power Consumption		
	1236 watts	1236 watts
Weight		
	468 lbs (175 kg)	468 lbs (175 kg)
BTU/Hr		
	4319	4319