



5085 Corporate Exchange Blvd. S.E.
Grand Rapids, MI 49512 • 616.656.2900

SS2201-J2
SS2201-J2G

SINGLE COMPARTMENT STAINLESS STEEL WARMING CABINET

Introduction

SKYTRON's single compartment warming cabinets are designed to maintain clinically ideal temperatures for sterile water, saline solutions and blankets in the care of patients in surgery, recovery, OB/GYN, ICU, ER, trauma areas and patient floors. The SKYTRON model SS2201-J2(G) can be countertop or under counter installed.

Technical Data

SKYTRON Warming Cabinets are manufactured from stainless steel with double wall insulation and polished to a #4 satin finish. The door is constructed of double pan stainless steel and insulated and polished to a #4 finish. Optional glass door is double pane tempered glass for low heat loss. Warming chamber designed to maintain temperatures through the use of digital controls. The SS2201-J2(G) has an interior compartment of 4.37 cubic feet.



Standard Features

- Silicone Door Gasket
- Cord and Plug (8 feet)
- Illuminated Power ON/OFF Switch
- Digital Thermostatic Control and Display
- Lock-out Capability (Temperature Setting)
- Reversible Door Swing
- Temperature Range 90°-160° F
- Fahrenheit Display (Celsius optional)
- Front Access Serviceability
- Leveling Feet

Options

- Solid door configuration SS2201-J2
- Glass door configuration SS2201-J2G
- Warmer cart F1-020-23
- Warmer stand F1-020-23-NC
- Caster Kit F2-012-58
- Celsius control

Exterior Dimensions: 30"W x 26"H x 26"D

Interior Dimension: 27"W x 15"H x 19"D

Electrical: 120V, 60Hz, Single Phase,
8-foot Power Cord with Plug

Power Consumption: 240 Watts

Amperage: Approximately 2.0 Amps

Weight: 179 lbs.

BTUs: 819/hr.

